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CRISPR/Cas9-Mediated Genome Editing in *C. elegans*: using Cas9 Protein + Purified sgRNAs)

-----Gang Wan's lab at School of Life Science, Sun Yat-sen University, China
(Version. 20250201)

Note: The following protocol was translated by ChatGPT from the Chinese version above.

1. sgRNA Design and Purification

1.1 sgRNA Design

Design sgRNAs targeting the desired region of the gene based on experimental goals:

- **(1) Knockout of the target gene**

Use SnapGene to search for GNGG motifs near the start and stop codons. Design 1 sgRNA for each end. Evaluate efficiency using CRISPOR (scores are for reference only).

Notes:

- GC content of the 20 bp upstream of NGG should ideally be 40–60%.
- Avoid poly-A/T/C/G stretches.
- If the first nucleotide is not G, add a G to the 5' end.
- Prefer guides targeting the sense strand, as antisense strand guides tend to be less efficient (ignore what the software says about strand).

- **(2) Deletion of a specific gene fragment**

Design 1 sgRNA upstream (5') and 1 downstream (3') of the region to be deleted.

- **(3) Site-specific base substitution**

Design an sgRNA within 10–20 bp of the target base. If it's too far, editing may occur at the sgRNA cut site but not at the desired base.

- **(4) Tag insertion into target gene**

Insert tags near the 3' end of the start codon or 5' end of the stop codon to allow for in-frame fusion. Design 1 sgRNA as close as possible (within 10–20 bp) to the insertion site.

1.2 sgRNA Amplification and In Vitro Transcription

1. Use pGW145-1 as the template. PCR amplify the sgRNA fragment (includes T7 promoter, sgRNA, and tracrRNA).

2. Perform in vitro transcription using T7 RNA Polymerase (NEB M0251L or Vazyme equivalent):

In a 0.2 mL tube on ice, prepare 20 µL reaction:

- 2 µL 10× RNAPol buffer
- 0.4 µL NTP mix (NEB N0466S)
- 280 ng sgRNA PCR product
- 0.5 µL RNase inhibitor (Takara 2313A)
- 0.25 µL 1 M DTT
- 2 µL T7 RNA polymerase
- RNase-free H₂O to 20 µL

Incubate at 37 °C for 4 h. Check product via 2% agarose gel.

3. Add:

- 10 µL DNase buffer
- 2 µL DNase I (Thermo Fisher EN0525)
- 68 µL RNase-free water

Incubate at 37 °C for 30 min to digest template DNA.

1.3 sgRNA Extraction

1. Add 500 µL TRIzol and 100 µL chloroform to sgRNA solution. Shake, incubate, and centrifuge at 18,000 × g, 4 °C for 10 min.

2. Transfer aqueous phase to new tube. Add 2.5× volume cold 100% ethanol. Precipitate overnight at –80 °C.

3. Centrifuge at $18,000 \times g$, 4 °C for 10 min. Wash pellet with 75% ethanol. Dry and dissolve in 12 μ L RNase-free water. Check RNA integrity via gel electrophoresis. Store at -80 °C.

2. Repair Template Construction and Purification

2.1 Template Construction

- **(1) Whole gene knockout:** No repair template needed.
- **(2) Deletion of specific fragment:** Clone ~100–500 bp homology arms flanking the target region into pGEM-7ZF(+) vector (EcoRI/HindIII sites). Sequence, then amplify as template.
- **(3) Base editing:** Include mutated base within homology arms.
- **(4) Tag insertion:** Insert tag coding sequence between homology arms at start or stop codon region.

2.2 Template Purification

Add 900 μ L cold E4 buffer (5mL Tris-HCl(1 M, pH8.5), 12.5 mL NaCl(5 M) for 50ml) + 660 μ L isopropanol to the recovered DNA. Mix, spin at $18,000 \times g$, 4 °C for 30 min. Wash with 70% ethanol, dry ≤ 10 min. Resuspend in 15 μ L RNase-free H₂O. Store at -20 °C.

3. Microinjection

3.1 Worm Preparation

Pick 20–50 L4-stage worms one day before injection. Let them grow overnight to adults.

3.2 Agarose Pad Preparation

Prepare 2.5% agarose. Drop 50 μ L onto a coverslip and overlay another slip in a cross (“+”) pattern. After 2 min, remove the top slip — pad is ready after 6 h at room temp.

3.3 Injection Mix Preparation

Thaw sgRNAs (for **unc-58** and target gene), Cas9 protein, repair templates, and RNase-free water on ice. Heat-denature dsDNA repair templates using a PCR program:

95 °C 2 min → 85 °C 10 s → 75 °C 10 s → 65 °C 10 s
→ 55 °C 1 min → 45 °C 30 s → 35 °C 10 s → 25 °C 10 s

Prepare injection mix (10 µL total):

Component	Amount
sgRNA (unc-58)	part of 1.8 µg total sgRNA
sgRNA1 (target gene)	-
sgRNA2 (target gene)	-
Cas9 protein	5 µg
RNase-free H ₂ O	to 10 µL

Incubate 15 min at 37 °C. Then add:

- AF-JA-76: 240 ng
- Repair template: 250 ng (if applicable)

Spin at 18,000 × g, 4 °C for 2 min. Transfer 7.5 µL of supernatant to fresh tube for injection.

3.4 Needle Preparation

Use WPI 1B100F-4 capillaries and pull with Narishige PC-100 to generate injection needles. Load ~2.5 µL injection mix. Let gravity fill the tip.

3.5 Needle Calibration

Use a heated capillary to make a glass fiber. Place under microscope with oil droplet. Use it to align and calibrate the injection needle with a reference. Set injection pressure to 0.5 and injection value to 0.6. Confirm injection works.

3.6 Worm Injection

Pick young adult worms (with early embryos), wash on OP50-free NGM plate. Mount 5 worms in oil on the agarose pad. Under 40× magnification, inject into the

middle of the germline twice per worm. Transfer injected worms into fresh NGM plate using M9 buffer.

4. Worm Screening and Genotyping

Label injected worms as **P0**. After 3 days, pick 96 F1 worms with *unc* phenotype, transfer to individual plates. After F2 appear, extract F1s and genotype using a worm lysis protocol (include adults + larvae to reduce false positives).

From F1 plates with desired heterozygous genotypes, pick 12 F2s (wild-type movement), grow for 2–3 days. Once F3s hatch, genotype F2s to identify homozygous mutants.

After identifying homozygotes, **outcross twice with wild-type (N2)** to eliminate off-target effects. Freeze validated lines.

Important: Always isolate at least 2 independent alleles and outcross before experiments or freezing.